

Transparent Film Island Dressing

Asguard Transparent Film Island Dressings provide a low adherent and breathable absorbent pad permeable to moisture. The waterproof film assists with unnecessary dressing changes and help protects from infection. The adhesive is skin-friendly and provides reliable fixation, yet is easy to remove, minimising trauma to the healing tissue.



Flexible

Ideal for use on areas difficult to dress areas (e.g. heel)




Shower Proof

Water resistant film allows showering with the dressing in-situ



Non-Adherent Pad

Comfortable and absorbent with minimal pain on removal



Suitable for:

- Post-operative wounds
- Minor cuts and abrasions
- Wounds with a low level of exudate

| REF | Dressing Size | Padding Size | Sterility | Pcs/Unit | Units/Ctn |
|---------|---------------|--------------|-----------|----------|-----------|
| ISC001 | 4cm x 5cm | 2cm x 2.5cm | Sterile | 50 | 10 |
| ISC002 | 6cm x 7cm | 3cm x 4cm | Sterile | 50 | 10 |
| ISC0610 | 6cm x 10cm | 2.5cm x 6cm | Sterile | 50 | 10 |
| ISC0910 | 9cm x 10cm | 4.5cm x 6cm | Sterile | 25 | 10 |
| ISC0915 | 9cm x 15cm | 4.5cm x 10cm | Sterile | 25 | 10 |
| ISC003 | 10cm x 12cm | 4.5cm x 6cm | Sterile | 25 | 10 |
| ISC028 | 10cm x 20cm | 5cm x 15cm | Sterile | 25 | 10 |
| ISC025 | 10cm x 25cm | 5cm x 20cm | Sterile | 25 | 10 |
| ISC026 | 10cm x 30cm | 5cm x 25cm | Sterile | 25 | 10 |
| ISC0935 | 10cm x 35cm | 5cm x 30cm | Sterile | 30 | 10 |
| ISC004 | 15cm x 20cm | 8cm x 10.5cm | Sterile | 25 | 10 |

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

