

seca 403 and seca 402

Carts for mobile support of seca baby scales

new



- Custom-made indentations in the surface hold seca baby scales firmly in place.*
- Ideal working height of 80 cm / 31.5 inches.
- Integrated drawer and shelf offer extra storage space so you always have everything within reach (seca 403).
- Optionally available paper roll holder helps you keep the work area clean (seca 402).

seca 403 and seca 402:

The perfect supplement to
seca baby scales.

The seca carts seca 403 and seca 402 are ideal supplements to the seca baby scales.* They were specially developed for safe and hygienic routine use in hospitals and medical practices. Customized indentations in the surface accommodate the scale's base and ensure that the scale does not slip out of place. Two practical brakes keep the cart stable so babies can be safely weighed. And you can push or pull the cart by its handle to move the baby scale quickly to the place where it's needed.



For work that's simple,
ergonomic and hygienic.

Mobile and safe use of baby scales is of prime importance in pediatrics. So now the seca baby scales* can be securely positioned on specially developed carts. The ergonomic design of the carts, with a work surface at 80 cm, ensures safe work processes. Hygienic standards are easier to maintain, thanks to the seamless, easy-to-clean surface.

The seca 403 offers practical
storage options.

The expanded variant of the seca cart has practical storage space that keeps documents and materials close at hand. Extra space is provided by a shelf and an integrated drawer.

Safe measuring and weighing
anywhere – with the seca 402.

In the basic model too the seca cart offers all the advantages of safe, mobile and hygienic work with seca baby scales. A practical and simply installed paper roll holder is optionally available.

The optional paper roll holder on the seca 402 makes work even easier and more hygienic.



Customized indentations keep the seca baby scale anchored on the cart.



Four bumpers on the wheels protect the cart and scale from collisions with walls and doors.

seca 403 and seca 402

Technical Data

seca 403:

- Dimensions (WxHxD): 715 x 800 x 435 mm / 28.2 x 31.5 x 17.1"
- Weight: 27.5 kg / 60.6 lbs
- Max. load: 100 kg / 220 lbs
- 4 bumpers
- 2 brakes
- Integrated drawer and shelf

seca 402:

- Dimensions (WxHxD): 715 x 800 x 435 mm / 28.2 x 31.5 x 17.1"
- Weight: 11.5 kg / 25.4 lbs
- Max. load: 100 kg / 220 lbs
- 4 bumpers
- 2 brakes
- Optional: paper roll holder

* Both carts offer an individual surface for following seca baby scales: seca 374, seca 376, seca 334, seca 336, seca 354, seca 382, seca 383, seca 384, seca 717, seca 728, seca 727 and seca 757.



☎ 1800 67 5432

📄 1800 67 5433

✉ sales@ecomед.com.au

🌐 www.ecomed.com.au