



# Welch Allyn Connex<sup>®</sup> ProBP<sup>™</sup> 3400 Digital Blood Pressure Device

Accurate, reliable, connected blood pressure readings  
in the palm of your hand



Take blood pressure accuracy to the next level.  
In fact, take it anywhere.



The Welch Allyn Connex ProBP 3400 puts consistent, accurate blood pressure readings at everyone's fingertips.



- Small, handheld design makes device easy to carry and capture readings on rounds or in triage settings
- SureBP® NIBP technology enables fast, accurate 15-second blood pressure readings, improving patient comfort
- Obtains accurate readings and minimizes error codes, even in the presence of many cases of motion or weak pulse
- USB-wired technology enables integration to select EHR/EMR systems, helping improve documentation accuracy and efficiency
- Compatible with Welch Allyn cuff sizes ranging from small child (size 8) to thigh (size 13)
- Lithium-ion battery technology delivers up to 100 readings on a single charge
- Configurable to match your workflow and specific patient needs. Upgradeable so that your initial investment is protected as your needs change

*Small, handheld design makes device easy to carry and capture readings on rounds or in triage settings*

Configured for your facility and your workflow.

Welch Allyn Connex ProBP 3400	Handheld	Wall Mount	Mobile Stand	Desk Mount
SureBP	34XFHT-B*	34XFWT-B*	34XFST-B*	4602-61 Desk Mount only. Requires Handheld ProBP (not included).
Standard BP	34XXHT-B*	34XXWT-B*	34XXST-B*	

\*Power-cord configurations: Replace "B" in part number with designation from chart below.

North America	B
Europe	2
UK	4
Australia (orange cord)	6
China	C
South Africa	7
Japan	J'

\*Mobile stand version only. Handheld version uses North American "B" configuration

## The perfect blood pressure solution for just about any care setting.



### Physician Office

With AAMI/BHS-AA-certified accuracy, SureBP® blood pressure technology that captures readings in just 15 seconds, and the ability to work with a broad range of cuff sizes<sup>1</sup>, the ProBP™ 3400 is a reliable alternative to manual sphygmomanometers. The ProBP 3400 is also an affordable approach to digitally capture blood pressure measurements and export readings to electronic health record systems.



### Hospital

Small enough to carry in your hand, the ProBP 3400 is a convenient alternative to manual sphygmomanometers or larger vital signs devices for capturing reliable blood pressure readings during rounds or in triage settings. A mobile stand featuring a compact footprint and accessory basket enhances portability with an added measure of security and the ability to store extra cuff sizes on board.



### Specialty Care

In situations when you need accurate blood pressure readings quickly and easily, look to ProBP 3400 first. Whether in long-term care, sports medicine, dental health, or other applications that call for reliable, repeatable blood pressure readings, the portable and accurate ProBP 3400 is the right choice for you.

*"...(SureBP® technology) offers great advantages in speed and patient comfort...clinicians can have added confidence in this new technology."\**

*-Dr. Bruce S. Alpert*



### Proven Accuracy and Reliability

Using the same ground-breaking SureBP® technology found in the Welch Allyn Spot Vital Signs LXi, the ProBP™ 3400 can capture blood pressure readings in just 15 seconds—nearly twice as fast as other electronic blood pressure devices.

- > Meets AAMI standard for blood pressure accuracy
- > Achieved "AA" rating by British Hypertension Society data analysis

\*Alpert, Bruce S., Clinical evaluation of the Welch Allyn SureBP® algorithm for automated blood pressure measurement. Blood Pressure Monitoring, Vol. 12 No. 4, August 2007.



## Connex® ProBP™ 3400 Digital Blood Pressure Device



### Connect to your electronic health record system.

The Connex ProBP 3400 comes connectivity-ready, allowing you to integrate it with the electronic health records system you have now or in the future. With hardwired USB connection, the Connex ProBP 3400 communicates with select EHR/EMR systems in both the physician's office and the hospital, helping reduce documentation errors and improve efficiency.



### The correct cuff means the correct reading.

Unlike many digital blood pressure devices, the Connex ProBP 3400 works with a very broad range of cuff sizes<sup>1</sup>, minimizing the potential for inaccurate readings on patients large or small. Plus you have all the advantages of unique Welch Allyn FlexiPort® cuffs that enable quick, snap-fit connections, disposable cuffs to reduce cross-contamination risk and facility-standardization benefits. Visit [www.welchallyn.com/flexiport](http://www.welchallyn.com/flexiport) for more information.

<sup>1</sup>The ProBP 3400 is not intended for neonatal use.



### Protect your investment

Remote troubleshooting. Software updates and upgrades. Extended warranties.

Go beyond the standard product warranty to take service and support to another level with Welch Allyn Partners in Care Services<sup>SM</sup>. Get faster turnaround times, free shipping, software upgrades, and more—everything you need to keep your device on site and in service longer.\*

- Comprehensive Partnership Program (1 yr) / S1-3400
- Comprehensive Partnership Program (2 yr) / S1-3400-2

\*Availability of specific features, coverage, hours and response times may vary by location. Other limitations may apply.

For more information about Connex® ProBP 3400 please contact your Welch Allyn representative or visit [www.welchallyn.com/connex](http://www.welchallyn.com/connex).

Welch Allyn Corporate Headquarters  
4341 State Street Road, P.O. Box 220  
Skaneateles Falls, NY 13153-0220 USA  
(p) 800.535.6663 (f) 315.685.3361

[WWW.WELCHALLYN.COM](http://WWW.WELCHALLYN.COM)

**WelchAllyn®**