



Low Linting



Anti-Static



Infection Control

Owear Surgical Hoods are constructed from a non-woven, spun bond polypropylene fabric. The ties are 42cm in length to allow for easy application and conform to the head for a comfortable fit.

TIE AT THE NECK

Adjustable tie for personalised fit

BREATHABLE

Comfortable for prolonged periods

COMPLETE COVERAGE

Pull over hood for overall protection



Suitable for:

Hospital and medical

Operating theatre

Delivery suites

ICU

Reorder	Feature	Colour	Pcs/Unit	Units/Ctn
HD001	Regular	Blue	250	4
HD002	Small	Blue	250	4